



Take Control

IQ/OQ/PQ Kits, Tools and Software
for Analytical Instrument Qualification

VKIT-LFM2 HPLC / UPLC flowmeter

VKIT-LFM2 HPLC / UPLC flowmeter

Specially designed for HPLC and UPLC, the V:Kit liquid flowmeter enables easy and accurate measurement of flow rates. The tripod mount and robust carry case make VKIT-LFM2, the perfect tool for field service engineers and laboratory personnel performing OQ/PQ and routine performance verification activities.

The V:Kit HPLC liquid flowmeter displays the current flow rate in mL/min to four decimal places. A USB port is provided for connection to a PC or serial printer, enabling compliance with the data integrity demands required in GxP environments.

Using the HPLC liquid flowmeter with the optional AutoFlow software or printer enables unattended operation: particularly useful when measuring the low flow rates associated with LCMS.

The flowmeter makes possible simpler diagnosis of common HPLC pump problems such as poorly functioning check valves or leaking pump seals.

LFM2 improves on the classic, well-proven VKIT-LFM design by adding faster response time, greater accuracy at low flow rates, a more compact design and improved connectivity via USB.

Technical specifications

- 0.05 to 25.00 mL/min flow range
- Flow accuracy $\leq 1.0\%$
- Flow precision $\leq 0.1\%$
- 0.1 $\mu\text{L}/\text{min}$ resolution
- Multi-point calibration with UKAS (UK equivalent of NIST) traceability calibration certificate
- PTFE, PEEK, FEP, EPDM & Acetal wetted parts
- 195 μL sensor volume: 5 readings each minute at 1.00 mL/min.
- 24V DC UK (3pin UK 240V), EU (2-pin 220V) or USA (2-pin 120V) (state when ordering)
- Backlit LCD display
- Weight: 1.8 kg
- Dimensions: 132 mm x 89 mm x 64 mm (HxWxD)
- Integrated tripod stand allows easy positioning on the bench
- Supplied with sturdy carrying case



Integrated tripod stand



Carry case

About V:Kit

V:Kit specialise in delivering solutions for Analytical Instrument Qualification. With a comprehensive range of measuring tools, software, protocols, standard operating procedures and certified standards, V:Kit products are used by OEMs, service providers and pharmaceutical laboratories worldwide.

Take control of your qualification needs with V:Kit.


GLP / Data Integrity Options

To assist with GLP and the increasing demands for data integrity by regulatory agencies such as FDA, EMA and PMDA, the liquid flowmeter can be optionally specified with software or a dedicated printer.

VKIT-PRT Printer

The flowmeter may also be specified with a printer option, to provide simple, effective traceability.

VKIT-PRT This is connected via the USB port and provides a time and date stamp of each flow reading. No PC software is required, and the printer is encoded with the flowmeter serial number.



AutoFlow

AutoFlow is an optional software module. It is available as standalone (flowrate only) or as part of the V:Kit 5 Analytical Instrument Qualification package.

Flow data can be read directly from the flowmeter into V:Kit AutoFlow software to produce an audit-friendly and traceable calibration report which evaluates the flow accuracy and precision of your HPLC and UPLC equipment.

V:Kit Take Control Software for Analytical Instrument Qualification		HPLC Flow Rate Certificate	
Certificate Number:	20210526-155043-12	Protocol:	HPLC Flow Pump Test
Customer Name:	V:Kit Ltd	System Location:	Agilent Park
Module Serial Number:	DE64301213	Module Model Number:	G1311A
Test Certificate:	123456	Test Serial Number:	2018-999
Calibration Date:	01-Jan-2021	Calibration Expiry:	01-Jan-2022
Qualifier:	V:Kit Support	Department:	Agilent 1100 HPLC SYS_03
Test Name:	HPLC Flow Pump Test	Test Module Type:	VKIT-LFM2
Test Type:	HPLC/UPLC Flowmeter	Test Date:	01-Jan-2021
Test Expiry:	01-Jan-2022		
Scope:	This test uses a traceably certified liquid flowmeter to qualify an HPLC pump for flow rate accuracy, precision and stability at various flow rates.		
Accuracy and Precision	Set Flow:	0.20 (mL/min)	0.50 (mL/min)
	Measured Flow:	0.2000	0.5000
	Reading 1:	0.2011	0.5000
	Reading 2:	0.2010	0.5000
	Reading 3:	0.2010	0.5000
	Reading 4:	0.2008	0.5005
	Reading 5:	0.2010	0.5004
	Mean Flow:	0.20096 mL/min	0.50015 mL/min
	SD:	0.0001028 mL/min	0.00028 mL/min
	Accuracy:	0.48 %	0.01 %
	Limit:	±2.00 %	±2.00 %
	Status:	PASS	PASS
	PSD:	0.05 %	0.25 %
	Limit:	±2.00 %	±2.00 %
	Status:	PASS	PASS
	Limit:	0.0000	0.0000
	Regression:	1.0000	0.99900
	Status:	PASS	PASS
Overall Test Status:	PASS		
Additional Documents Related to this Test:	Flowmeter Certificate of Calibration		
Role:	Name:	Signature:	Date:
Qualifier:	V:Kit Support	Electronically signed by V:Kit Support	26-May-2021

Typical Flowrate Report



VKIT-LFM2 HPLC / UPLC flowmeter

More about V:Kit 5.1 Software - Qualify your laboratory systems with a 21 CFR Part 11 compliant package

- Generate comprehensive qualification documents with full traceability.
- Create, manage, and distribute secure, version controlled executable qualification protocols directly linked to Standard Operating Procedures and Certificates of Analysis.
- Manage your Analytical Equipment database.
- Share qualification results and protocols with your team. User account management with password expiry and security roles.
- Electronic Signatures, Audit Trail, Deviation Reporting.

Available for HPLC/UPLC, GC/GCMS, LCMS and MS/MS, Dissolution and UV-Vis.

Ordering codes

VKIT-LFM2	V:Kit HPLC/UPLC flowmeter, with calibration certificate. (Specify your required number of calibration setpoints, 1,2,3,4 or 5. Power supply requirements UK, EU, US or other)
V5-2002	V:Kit HPLC AutoFlow Software Starter Bundle
V5-PRT	Printer for VKIT-LFM/LFM2 HPLC flowmeter

Ordering information

Order online at www.v-kit.com. We accept payment with VISA, MasterCard and Bank Transfer in EUR, USD and GBP currencies.

Or send your Purchase Order to sales@v-kit.com.